

Shelf Option for SLS PCS

Convenient, Durable, Reliable

Designed to fit into your existing workflow, Codonics Safe Label System (SLS) optional shelf puts the SLS Point of Care Station (PCS) at your clinician's fingertips and often where space is very limited. The small, sturdy footprint enables SLS PCS to be easily sited into your operating room or anywhere medications are prepared.

To maximize convenience and user preferences, the shelf's height is adjustable* and it can be conveniently mounted on the left or right side* of your cart. The shelf includes a locking drawer that can be used for additional SLS PCS supplies to be readily available at the point of care. The shelf includes mounting hardware designed to securely hold SLS PCS in place. Ask your Codonics sales representative for details.

Codonics shelf is compatible with the following manufacturers' carts:

- Armstrong Medical Industries Inc.®
- Capsa Healthcare®
- Waterloo Healthcare® Medical Carts
- Blue Bell Medical®
- Metro®
- Omnicell®



KEY FEATURES

Advantages of the SLS PCS Shelf

- Enables placement wherever medications are prepared
- Added storage in locking drawer holds ink and labels as well as other supplies**
- Cables are easily managed and secured
- Integrated power supply bay
- Chassis resistant to disinfectants and hospital-grade cleaners
- Configurable for left or right side mounting

Drawer Dimensions: 10.7" (27.1 cm) W,
12.1" (30.7 cm) D, 2.7" (6.8 cm) H

Shelf Dimensions: 11.1" (28.2 cm) W,
15.7" (39.9 cm) D, 3.8" (9.7 cm) H***

Total Weight: 13.5 lbs. (6.1kg)
(varies depending on mounting option)

* Dependent on the cart manufacturer's recommendations

** Some PCS media cannot be accommodated

*** 11.1" (28.2 cm) H, including support bracket



We bring the future into focus

17991 Englewood Drive
Middleburg Heights, OH 44130 USA
+1.440.243.1198
+1.440.243.1334 Fax
info@codonics.com
www.codonics.com

Call Codonics today at **800.444.1198** or visit www.codonics.com for more information.

All registered and unregistered trademarks are the property of their respective owners. Specifications subject to change without notice. Patents: www.codonics.com/ip/patents Printed in the U.S.A. Copyright © 2017 Codonics Inc. 12/2017



Sponsoring Donor

Corporate Member

900-603-003.02