

Your Image & Data DISTRIBUTION HUB



NEW FEATURE

Copy to USB allows users to record patient studies, reports and viewers onto USB drives or cards.



Scheduled Archive automatically records all studies to disc, including an incremental log file, not archived since the last session.



Built-in web server allows users to remotely access and control all functions of Virtua from a web browser.



User Interface is available in multiple languages.



Automatically record DICOM exams from modalities without tying up workstations.



Attach Non-DICOM files such as Word, PDF, MPG, JPG.



Encrypt patient data using AES 256 for maximum security.



Record over 62 CDs or 31 DVDs an hour with brilliant, full-color disc labels.

----- NETWORK -----



NEW FEATURE

Kiosk Mode query/retrieve patient studies from PACS with configurable restrictions to protect patient data.

Send studies from PACS to Virtua, as well as reports from RIS in DSR or HL7.



Direct to Disc records files (PDF, PPT, AVI, MPG) to disc using standard Windows file copy as an alternative to DICOM store including control over disc content and labeling.