


Codonics Clarity[®] Viewer

Codonics Clarity Viewer features simple image navigation and selection, an intuitive user interface, quick viewer launch and rapid image loading.



Codonics Virtua[®] Medical Disc Publisher is the new standard for speed, efficiency and ease of use in automatic disc recorders. Virtua allows you to include the industry's leading viewers while its brilliant, full color disc labels increase your brand awareness.

For CT, MR,
Ultrasound, PET**

- General DICOM image viewing includes full screen or mosaic mode with up to 16 x 16 images
- Window leveling, color LUTs, magnification, and measurement tools
- Windows printing via icon-based intuitive user interface
- English, French and Japanese user interface. Unicode display of characters including English, French, Flemish, Italian, German, Swiss, Malaysian*, Ukrainian, Spanish, Portuguese, and Danish
- Slice navigation, zoom, synchronization, and pan

* patient name only

**not intended for primary interpretation

Call Codonics today at **+1-440-243-1198** or visit www.codonics.com for more information.

Not available on
-A Virtua models.



Corporate Headquarters:
17991 Englewood Drive
Middleburg Heights, OH 44130 USA
(440) 243-1198
(440) 243-1334 Fax
Email info@codonics.com
www.codonics.com