


# Siemens syngo® fastView from Codonics

Siemens syngo fastView is a stand-alone viewer for DICOM images provided on DICOM exchange media, and can be used on any Windows PC. Its operating concept is based on the easy-to-use syngo philosophy, learn one – know all.



Codonics Virtua® Medical Disc Publisher is the new standard for speed, efficiency and ease of use in automatic disc recorders. Virtua allows you to include the industry's leading viewers while its brilliant, full color disc labels increase your brand awareness.

For PACS and  
general radiology\*

- Layouts: 1:1, 4:1, 16:1, 1:2 vertical, 9:1 stripe and stack modes
- Special ultrasound movie mode for cardiac load ultrasound examinations
- Special 3D MPR mode for displaying volume datasets in three orthogonal views
- Special fusion mode for displaying pre-registered, fusion datasets side by side
- Save images in JPEG, BMP or AVI formats
- Multi-language support (English, German, French, Spanish, Chinese, and Japanese)

\* not intended for primary interpretation

Call Codonics today at **+1-440-243-1198** or visit [www.codonics.com](http://www.codonics.com) for more information.

Not available on  
-A Virtua models.



Corporate Headquarters:  
17991 Englewood Drive  
Middleburg Heights, OH 44130 USA  
(440) 243-1198  
(440) 243-1334 Fax  
Email [info@codonics.com](mailto:info@codonics.com)  
[www.codonics.com](http://www.codonics.com)